EN

DE

IT

FR

ES

MiroFix[®] Inject

Instructions for use
Gebrauchsanleitung
Istruzioni per l'uso
Mode d'emploi
Instrucciones de uso

Manufactured by:

aap Biomaterials GmbH & Co. KG
Lagerstraße 11-15
64807 Dieburg
Germany

Distributed by:

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File: 191/01/2007

Instructions for use

Radiopaque plastic for use in percutaneous vertebral augmentation

Purpose and properties

MiroFix® Inject is a radiopaque, injectable plastic for use in percutaneous vertebral augmentation during vertebroplasty or kyphoplasty. Mixing the two separate sterile components produces a ductile, injectable bone cement which, after hardening in situ, is able to stabilise neoplastic vertebral lesions and vertebral compression fractures. MiroFix® Inject cement powder also contains insoluble zirconium dioxide as an X-ray contrast medium.

Composition and pack size

1 sachet (24.00 g of powder) contains: poly(methyl methacrylate) 10.95 g; poly(methyl acrylate/methyl methacrylate) 1.75 g; zirconium dioxide 10.80 g; benzoyl peroxide 0.50 g.

1 ampoule (10 ml of liquid) contains: methyl methacrylate 9.93 ml, N,N-dimethyl-p-toluidine 0.07 ml, stabilised with 60 ppm hydroquinone.

Uses

During percutaneous vertebral augmentation, MiroFix® Inject is indicated for the fixation and filling of symptomatic vertebral compression fractures, vertebral tumours and vertebral haemangioma.

Contraindications

The use of MiroFix® Inject is contraindicated in the event of

- acute or not completely healed local infections and active systemic infections
- treatment-refractory coagulopathy or during anticoagulation
- pulmonary insufficiency or severe heart and lung disorders
- vertebra plana or if safe percutaneous access to the vertebra is not guaranteed
- unstable vertebral fractures with involvement of the posterior margin
- previous damage to the pedicle wall (transpedicular access)
- fractures or neoplasia with narrowing of the spinal canal (more than 20%) and myelopathy
- retropulsing vertebral fragments with myelopathy
- acute traumatic fractures of non-osteoporotic vertebrae
- a satisfactory response to conservative treatment
- asymptomatic stable vertebral fractures
- known hypersensitivity or allergy to one of the components of MiroFix® Inject

MiroFix® Inject is not suitable for use in arthroplasty.

Use during pregnancy and lactation

MiroFix® Inject must not be used during pregnancy or lactation.

Warnings relating to the age of the patients to be treated

There are no adequate studies of the use of acrylate cements in children. As the possibility of acrylate cements having adverse effects on bone growth cannot be ruled out, MiroFix® Inject is not recommended for use in children and people who are still growing.

In the absence of long-term data on the use of acrylate cements in vertebral augmentation, MiroFix® Inject is not recommended for use in younger patients (< 55 years old).

Side effects

Immediately after implantation of acrylate bone cements, the rise in pressure in the medullary canal can lead to a temporary drop in blood pressure. In rare cases, pulmonary embolism and myocardial infarction are also observed. These cardiovascular and respiratory side effects, which are known as **implantation syndrome**, result mainly from bone marrow constituents infiltrating the venous system.

The following additional undesirable effects have occurred after using acrylate cements: a temporary drop in blood pressure, elevated serum levels of gamma-glutamyl transferase (gamma-GT) up to 10 days after surgery, thrombophlebitis, heterotopic ossification, cardiovascular reactions such as temporary heart rhythm disorders, short-term cardiac conduction irregularities, arrhythmia, myocardial infarction and cardiac arrest, hypoxaemia, bronchospasm, pulmonary embolism, apoplexy.

There have been rare reports of hypotension with anaphylaxis including anaphylactic shock associated with cardiac arrest and sudden death.

Possible side effects with a direct link to vertebral augmentation with PMMA bone cements include damage to paravertebral structures caused by leaking cement. Extravasation of cement can cause: compression of the spinal cord or spinal nerve roots with paralysis as well as paraesthesia and loss of sensation, radiculopathy, nerve root pain, sciatica, intercostal neuralgia, neuritis, dysphagia, muscle or other soft tissue pain, abdominal damage/ileus, displacement of vessels with the risk of embolism, pulmonary embolism and cerebral ischaemia. Other possible complications of percutaneous vertebral augmentation are fever, deep and superficial infections, pneumonia, lung infection, pedicle fractures, rib fractures, the collapse of an adjacent osteoporotic vertebra, new vertebral fractures, bleeding and haematoma, procedural pain and burns to the skin caused by imaging procedures.

Monomer vapours can irritate the respiratory tract and eyes and may damage organs.

Interaction

No known interaction.

Information for use

MiroFix® Inject should only be used by doctors who are fully proficient in the techniques required for percutaneous vertebroplasty and kyphoplasty.

In addition, before using MiroFix® Inject for the first time, surgeons should familiarise themselves with the mixing and application process. A trial run of the mixing process is recommended. If special mixing and application techniques are to be used, the surgeon must read the relevant instructions first.

The protective outer packaging (aluminium/PE sachet) and blister pack containing the ampoule should be removed from the carton in the non-sterile area. After removal from the protective outer packaging, the polyethylene sachet ("peel-off sachet"), the outside of which is not sterile, and the blister pack containing the ampoule are opened one after the other in sterile conditions by a circulating nurse and passed to a member of the surgical team in the sterile area. When the vertebroplasty or kyphoplasty apparatus is in the correct position, MiroFix® Inject can be injected into the vertebra using either a conventional two-part (!) syringe or another suitable application system (details of how to do this can be found in the instructions for the system used). The ampoule is opened by breaking its neck. The inner sachet containing the cement powder is cut open using sterile scissors.

Dosage

The amount of MiroFix® Inject required for injection depends on the anatomical conditions of the patient. However, at least one complete unit (the contents of one sachet and one ampoule) must always be mixed. Each dose should be mixed separately. It is advisable to keep additional packs of MiroFix® Inject ready as a precaution.

Tools required for mixing MiroFix® Inject

Sterile working area, sterile porcelain bowls, stainless steel bowls or plastic bowls that are suitable for monomers, sterile mixing spoons or spatulas made of porcelain or high-grade steel, or a sterile mixing system for bone cements.

Manual mixing

To make the mixture, empty the entire contents of an ampoule into a suitable sterile, inert mixing vessel. Then add the entire

contents of the pack of powder to the liquid and, using a suitable spatula, carefully mix the components for about 30 seconds until the powder has completely absorbed the liquid. The cement is ready for application when it has reached the desired viscosity. The handling time ends when the viscosity is so high that it is no longer possible to inject the cement. If the cement continues to be used after this point, there is a risk that the treatment will fail.

Using a mixing system

MiroFix® Inject can also be prepared in a mixing system. It does not need to be pre-chilled. The cement is designed for handling between 18°C and 23°C as standard and should be mixed for about 30 seconds. The contents of an ampoule are placed in the mixing vessel and all of the cement powder added. Please refer to the instructions supplied by the manufacturer of the mixing system for the rest of the procedure.

Application in a conventional two-part syringe

Once it has been mixed, the cement can be transferred to a two-part syringe and, after a given period (see Figure 1), injected under slight pressure. However, the procedure must be monitored closely by the surgeon. Please note that the injection time decreases as the size of the syringe increases, the thickness of the cannula decreases and the temperature increases. Conversely, the injection time increases as the size of the syringe decreases, the thickness of the cannula increases and the temperature decreases.

Caution: standard three-part syringes must not be used, as the rubber plunger that seals the barrel is unstable against the monomer component in MiroFix® Inject.

Application in special vertebroplasty or kyphoplasty apparatus

A few minutes after mixing, the cement can be transferred to special vertebroplasty or kyphoplasty apparatus and injected under slight pressure. In this case, the surgeon must read and follow the instructions supplied by the manufacturer carefully. Please note that the injection time for MiroFix® Inject also changes in relation to pressure (syringe size, cannula thickness) and temperature if special apparatus is used.

All application procedures:

After injection of MiroFix® Inject into the vertebra, the patient must maintain the same body position until the cement has hardened (see hardening times in Figure 1). This is the only way to achieve the best possible results.

Important general information

- Please note the temperature-time graph (see Figure 1).
 - Handling time and polymerisation are heavily dependent on the temperature of the components and environment. Hardening time is reduced by higher temperatures and prolonged by lower temperatures. The viscosity increases with the progression of polymerisation, i.e. the duration of the handling phase.
 - The addition of any other powders or liquids can reduce solidity and change the handling characteristics.
 - The powder and liquid components have been carefully designed to complement each other. The entire contents of the sachet and ampoule must always be mixed.
- USING ONLY PART OF THE COMPONENTS IS NOT PERMITTED.**
- Large air pockets should be avoided when filling the application apparatus.
 - Injection of MiroFix® Inject must be monitored using imaging procedures (real time).

Precautions and warnings

Percutaneous vertebroplasty and kyphoplasty should only be performed in medical facilities that can also perform open decompression in an emergency.

The long-term effects of acrylate cements on the vertebral column are unknown. Therefore, the use of MiroFix® Inject in younger patients (< 55 years old) is definitely not recommended.

There are no adequate studies of the effect of acrylate cements on human fertility.

Appropriate imaging procedures should be used to check that the vertebroplasty or kyphoplasty apparatus is in the correct position, that there is no damage to adjacent structures and that the injected MiroFix® Inject is in the right place. In addition, vertebroplasty and kyphoplasty should always be monitored using imaging procedures in real time. In particularly risky cases, the use of computed tomography as a further precaution should be considered.

During and immediately after injection of the bone cement/implantation, blood pressure, pulse and respiration must be monitored carefully and appropriate measures taken if significant changes occur. If the patient develops pulmonary or cardiovascular symptoms, appropriate monitoring of blood loss is required. In the event of acute respiratory failure, anaesthesiological measures should be initiated immediately.

In some circumstances, MiroFix® Inject can escape from the area of application in the vertebra or into blood vessels. The risk is increased if the bone cement is injected too quickly or under too much pressure, or if the vertebra is overfilled. If MiroFix® Inject escapes from the vertebra, the injection must be terminated immediately.

MiroFix® Inject must not be injected before the times specified in Figure 1, as the lower viscosity increases the risk of the cement escaping from the vertebra.

If the posterior margin of the vertebra is fractured, it is necessary to proceed with the utmost caution in order to prevent cement from escaping into the intervertebral foramina and vertebral foramen.

Based on current knowledge, filling the vertebra completely is not a prerequisite for effective pain relief. In fact, injecting too much cement increases the risk of uncontrolled leakage of cement from the area of application.

If the viscosity of MiroFix® Inject during the injection is too high, the apparatus and treatment can fail.

Should it be necessary to administer a second injection into the same vertebra, the apparatus used in the first application should be left in place in order to prevent the bone cement from escaping through the first puncture hole during the second injection.

Polymerisation (hardening) of MiroFix® Inject is an exothermic reaction. The heat given off during this reaction may damage bone or other tissue in the region of the implant.

The training and experience of the surgeon are very important when handling acrylic bone cements in vertebroplasty and kyphoplasty. Instructions for handling and mixing the cement must be followed carefully. Before using MiroFix® Inject, the surgeon must be thoroughly familiar with its properties and handling characteristics. As the handling and setting characteristics of MiroFix® Inject depend on the temperature and mixing technique, they are best determined by the surgeon on the basis of actual experience. For this reason, the surgeon is strongly recommended to carry out a trial run of the entire mixing, handling and setting process before performing a surgical procedure with MiroFix® Inject.

Methyl methacrylate is a slightly volatile, flammable liquid. The vapours produced during the mixing process can irritate the respiratory tract and eyes and cause general malaise and headaches. Such symptoms can be reduced with adequate ventilation or by using closed mixing systems. The monomer (methyl methacrylate) is lipid-soluble. Direct skin contact with the liquid monomer should be avoided as far as possible, as allergic reactions (contact dermatitis) cannot be ruled out. It is therefore advisable to wear an additional pair of polyethylene (PE) gloves under normal surgical gloves when handling the cement. The following materials have also proved to be suitable for protective gloves: PVP (polyethylene, ethylene vinyl alcohol, polyethylene) and viton butyl.

Contact lenses should be protected from the escaping monomer vapours.

Storage

MiroFix® Inject should be stored away from direct sunlight and below 25°C (77°F). Do not use MiroFix® Inject after the expiry date.

Shelf life/sterility

The expiry date is printed on the outer carton, protective aluminium packaging and inner sachet. MiroFix® Inject must not be used after this date. The contents of unused open or damaged packs must not be resterilised and should therefore be discarded. The polymer powder must not be used if it shows yellow discolouration. The liquid monomer, the ampoule itself, the inside of the blister, the cement powder and the inner PE/paper sachet are all sterile. The product must not be used if the packaging is damaged.

Safe disposal

Mixed cement should set before it is discarded with hospital waste.
For separate disposal of powder and liquid, please ask local waste disposal authorities.

Information

For further information, please contact your supplier or the manufacturer directly.

Symbol definitions

STERILE IEO

Sterilised with ethylene
oxide

Not for reuse

STERILE A

Sterile filling

Caution! See instructions for use.

LOT

Batch designation

Use by

Inflammable


Maximum storage temperature

Protect from sunlight

Figure 1 Temperature–time graph for handling MiroFix® Inject

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