

MiroFix[®] Inject.



besser

Schmerzreduktion ohne Risiko einer Hitze-Nekrose

better

pain reduction without risk of a heat necrosis

mejor

menos dolores sin el riesgo de una necrosis del calor

migliore

meno dolori senza il rischio di necrosi del calore

meilleur

réduction de douleur sans risque d'une nécrose de la chaleur



Exzellent.

exzellent



kurze Anmischphase, optimale Viskosität, hohe Kompressionsfestigkeit

excellent

short mixing phase, optimal viscosity, high compression firmness

excelente

fácil para mezclarse, viscosidad óptima, alta firmeza



© PD Dr. K. Wilhelm, Bonn

eccellente

facile mescolarsi, viscosità ottimale, alta fermezza

excellent

phase de mélange courte, viscosité optimale, fermeté élevée de compactage

New

optimale Kontrolle

50% Zirkoniumdioxid für optimale Röntgenkontrolle

optimal check


50% zirconium dioxide for optimal x-ray check

control óptimo

50% del dióxido del circonio para una visibilidad óptima en rayos-X

controllo ottimale

50% di diossido dello zirconio per un controllo ottimale del raggio-X



© PD Dr. K. Wilhelm, Bonn

contrôle optimale

50% du dioxyde du zirconium pour un contrôle de rayon-x optimale

MiroFix[®] Inject

Verkaufseinheit/Quantity Box. 1 Stück/1 pcs

Artikel-Nr./Product Code. 07010310

PZN. 6059477

Vertreiber/Distributor.

MiroMed[®] Pfeleiderer GmbH, Frankfurt, Germany

Spezifikationen/Specifications.

50% Zirkoniumdioxid/Zirconium dioxide

Kompressionsfestigkeit/Compression firmness 1220 bar

Anmischphase/Mixing phase 30s

(23°C/73,4°F Raumtemperatur/Air temperature)


Aushärtetemperatur/Hardening temperature 70°C/158°F

Besser.

Optimal.


Pfeleiderer Medical Produkte GmbH

MiroFix[®] Inject.



MiroCan[®] Diagnostic

Abgestimmt auf MiroCan[®] Diagnostic 2 in 1 Vertebroplastiekanülen und MiroJect[®] Applikationssystem

Optimized for MiroCan[®] Diagnostic 2 in 1 Vertebroplasty cannulas and MiroJect[®] application system

Optimizado para MiroCan[®] Diagnostic 2 en 1 cánulas para la vertebroplastía y para MiroJect[®]

Ottimizzato per MiroCan[®] Diagnostic 2 in 1 cannulas per il vertebroplastía e MiroJect[®]

Optimisé pour MiroCan[®] Diagnostic 2 en 1 de canules de la vertébroplastie et le système d'application de MiroJect[®]

Zubehör.

MiroMed[®] Pfeleiderer GmbH.

Der Spezialist für radiologische Produkte.

The specialist for radiological products.

El especialista en productos radiológicos.

Lo specialista per prodotti radiologici.

Le spécialiste pour les produits radiologiques.

MiroMed[®] Pfeleiderer GmbH

Address. Lindleystr. 12
60314 Frankfurt | Germany

Internet. www.mirofix.eu

Internet. www.miromed.de

E-Mail. mail@miromed.de

Fon. +49 69 489878-0

Fax. +49 69 489878-70



MiroFix[®] Inject.

**Der moderne
Knochenzement für die
Vertebroplastie**

**The modern bone
cement for vertebro-
plasty**

**El cemento moderno
del hueso para la ver-
tebroplastía**

**Il cemento ortopedico
moderno per il verteb-
roplasti**

**Le ciment orthopédique
moderne pour la ver-
tébroplastie**



**Dauerhafte Stabilisie-
rung von Wirbelkörpern**

**Lasting stabilisation of
vertebral bodies**

**Estabilización duradera
de cuerpos vertebrales**

**Stabilizzazione durevole
degli enti vertebrali**

**Stabilisation durable
des colonnes vertébrales**

